

Surearly Midstream

Extra-sensitive pregnancy test with **10mIU/mL** sensitivity for early diagnosis

Surearly EARLY SIGN PREGNANCY TEST

New Innovation in Digital Healthcare
SUREARLY

Surearly Early Sign

Pregnancy Test

EASY-TO-READ

- Result in 3-5 minutes
- Detect hCG hormone the lowest levels

*Over 99% ACCURATE

Extra Sensitive
10mIU/mL hCG


1 test

* At detecting typical pregnancy hormone levels.

Highly sensitive detection of hCG at 10mIU/mL concentration. Surearly Early Sign Pregnancy Test shows the positive result much sooner.

Key Features

- Extra-sensitive (sensitivity level of 10mIU/mL)
- Early detection of pregnancy
- Easy to use midstream type
- Rapid result in 3-5 minutes
- hCG hormone detection in urine



Surearly Early Sign Pregnancy Test detects hCG at concentrations as low as 10mIU/mL while most home urine pregnancy tests detect hCG at 25-50mIU/mL. As Surearly Early Sign Pregnancy Test has a very high sensitivity, it can tell your positive result much sooner.

Single-use (1 test/box)

- Contents: One Test Stick
- Quality Certificate : CE0123 (certified by TÜV SÜD)



SUGENTECH

Global In-vitro diagnostic
Total platform leader



Lab Test

Multiplex Immunoblot Assay



POC Test

Point-Of-Care Test
for Multiplex Parameters



Digital Healthcare

Mobile Integrated
Personal Healthcare



Contact Us

Head office / R&D Center
Sejong Campus PIUM, 75,
Jiphyeondong-ro, Sejong-si,
Republic of Korea

Manufacturing Site
Jeongjungyeonje-ro,
Osong-eup, Heungdeok-gu,
Cheongju-si, Chungcheongbuk-do,
Republic of Korea

Sales Office
5F, 12, Yeongdong-daero 96-gil,
Gangnam-gu, Seoul,
Republic of Korea

Asia Pacific: +82-44-850-9760
Other Regions: +82-2-2014-5621
sales@sugentech.com

WWW.SUGENTECH.COM/EN

SUREARLY



SUGENTECH

SUGENTECH

Surearly SMART

Digital Healthcare Platform

Surearly SMART


Mobile healthcare platform for Female Hormone Diagnosis!
Rapid detection of various hormones such as for pregnancy(hCG), ovulation(LH), menopause(FSH) and others(Estrogen, Progesterone) by using the unique semi-quantitative analysis method.

Key Features

• Bluetooth Connection

- The Smartphone Application will provide calendar, result tracking, hormone pattern analysis and related advice, etc.

• Semi-quantitative Analysis : Hormone Level Detection




• More Cost-efficient

- Durable, Reusable Digital Device by battery replacement
- Users can keep the digital device and buy refill test sticks.

Why testing Hormone?

Hormones are an essential factor for life and health. As aging and stress can cause hormone imbalance, consistent analysis and care are important.



Surearly Digital

Korea's First & Best Digital Test

Surearly DIGITAL PREGNANCY TEST

ACCURATE and SIMPLE!
Test result is displayed on LCD with digital word and symbol.

Key Features

- Rapid, Easy-to-Read DIGITAL results in about 3 minutes
- Over 99% accurate
- User's procedural error detection
- hCG hormone detection in urine
- Sensitivity level of 25mIU/mL



Multiple-use (7 tests/box)


- More cost-effective
- Contents : One Digital Test Reader with 7 Test Sticks
- Quality Certificate : CE0123 (certified by TÜV SÜD)



Surearly SMART Pregnancy

hCG Hormone Detection for Pregnancy

- hCG hormone level detection
- Test on the day your period is due.



Surearly SMART Ovulation

LH & Progesterone Detection for Ovulation For More Fertile Days, Predict(High, peak) & Confirm


- More fertile days by dual-hormone detection
- Test at the same time each day and normally the morning.



Surearly SMART Menopause

FSH Hormone Detection for Menopausal Transition & Menopause

- Test 2 weeks after the initial test
- A first morning urine is recommended.




Surearly SMART Hormone DUO

Dual-Hormone Detection for Estrogen & Progesterone

- Health & Beauty Care by understanding the key female hormone pattern (i.e. skin care, emotional care, exercise, etc. during her cycle)
- Test everyday with your first morning urine.

Surearly SMART Mobile App

- Monitors personal hormone balance
- Manages personal symptom and condition
- Structurally manages life pattern and interest




Surearly DIGITAL OVULATION TEST

EASY-TO-READ and ACCURATE!
Test result is displayed on LCD with digital word and symbol.
[YES means ovulation will soon occur within the next 24~36 hours.]

Key Features

- Rapid, Easy-to-Read DIGITAL results in 3 minutes
- Over 99% accurate
- User's procedural error detection
- LH hormone detection in urine
- Sensitivity level of 25mIU/mL



Multiple-use (7 or 20 tests/box)

- Contents : One Digital Test Reader with 7 or 20 Test Sticks
- Quality Certificate : CE0123 (certified by TÜV SÜD), FDA Registered